

# MADISON EVO

TOTAL KNEE PROSTHESIS



S U R G I C A L   T E C H N I Q U E

**serf**  
KNEE

# INDICATIONS

**The Madison total knee prosthesis is intended to replace the knee joint and the indications are:**

- Osteoarthritis (primary and secondary);
- Revision of a unicompartmental knee prosthesis.



**Madison instrumentation is designed to adapt to a traditional surgical approach.**

The surgeon is guided by bony landmarks to ensure optimal implant positioning.

Different bone anchoring techniques (cemented or cementless) can be used for femoral and tibial components.

**CONTENTS**

|                             |    |                            |    |
|-----------------------------|----|----------------------------|----|
| IMPLANTS .....              | 4  | TRIALS .....               | 29 |
| DESIGN .....                | 5  | TIBIA .....                | 31 |
| MECHANICAL TESTING .....    | 6  | EXTENSION ROD OPTION ..... | 32 |
| SURGICAL TECHNIQUE .....    | 7  | PATELLA .....              | 33 |
| INCISION AND EXPOSURE ..... | 8  | IMPLANTS .....             | 34 |
| TIBIAL PREPARATION .....    | 9  | IMPLANT DIMENSIONS .....   | 35 |
| FEMORAL PREPARATION .....   | 18 | IMPLANT SIZES .....        | 36 |
| SPACERS OPTION .....        | 20 | INSTRUMENTATION .....      | 38 |
| FEMUR PREPARATION .....     | 21 |                            |    |



# IMPLANTS



## Femoral components

- 7 **MADISON EVO** sizes and 1 **MADISON** size
- Cobalt chrome
- Posterior stabilised (PS) or cruciate retaining (CR)
- HA-coated cementless or cemented
- Anatomical trochlea

## Tibial inserts

- 8 sizes (1 to 8)
- UHMWPE high density polyethylene
- 4° built-in posterior slope
- Thicknesses of 10 to 20 mm (2 mm inc.) for posterior stabilised (PS) and Ultracongruent (UC)\*\*

## Tibial baseplates

- 8 sizes (1 to 8)
- Cobalt chrome
- Symmetrical
- Mobile or fixed bearing
- HA-coated cementless or cemented

## Patella components

- 4 diameters: 30, 33, 36 and 39 mm
- 2 thicknesses: 8 and 10 mm
- UHMWPE high density polyethylene
- Resurfacing and cemented with 3 pegs

## Optional Tibial extension rods

- Titanium alloy
- Diameter 11, L 35 & 55 mm
- Diameter 11, L 95 mm\*
- Diameters 9 & 13, L 35, 55 & 95 mm\*



\* Please check product availability with your local representative.

\*\* Some thicknesses are optional. For more details please contact your local representative.

# DESIGN

## **Madison implants are designed to preserve the bone stock.**

- 8 mm thickness for the femoral component (distal and posterior cuts)
- Trochlea with smooth anterior internal surface
- Minimal resection for the PS cage
- V-shaped fluted keel to maximise the contact area, without macrostructure and easy to remove in the event of revision surgery

## **Madison implants are designed to ensure long-term results.**

- Large contact areas
- Single ML radius
- Severe mechanical testing

## **Madison instrumentation is designed to be precise.**

- AP and ML size and external femoral rotation on the same instrument
- Real external femoral rotation
- Whiteside's line angle finder
- Unique design for the validated tibial extramedullary option with navigation
- Micrometric adjustment on the tibial cutting guide available for extramedullary or intramedullary alignment options

## **MADISON instrumentation is designed for a standard procedure.**

- Optimised dimensions for each instrument:
  - Medialised distal femoral cutting guide
  - Medialised proximal tibial cutting guide with patellar tendon protection



# MECHANICAL TESTING

## Wear testing according to ISO 14243-1 and 14243-2

4 prostheses (femoral component, UHMWPE insert and tibial baseplate) were installed on simulators reproducing gait cycles.

Tests were carried out at a frequency of 1 Hz over 5 million cycles.

- Results: all 4 prostheses passed the test without delamination or UHMWPE insert damage or fracture.

## Dynamic tests conducted on tibial baseplates according to standards ISO 14879-1 and ISO 21536

5 tibial baseplates were tested in flexion at 10 Hz over 10 million cycles, under a 900N compression load.

- Results: all 5 tibial baseplates passed the test with no distortion, breakage or damage.

## Shear tests conducted on six femoral components

These tests were conducted on **MADISON** femoral components as there is no reinforcement under the trochlea, thus avoiding an additional cutting stage during surgery and preserving the bone stock.

During the tests, one condyle of each femoral component was fixed and the load was applied to the opposite condyle.

- Results: 5 femoral components broke under a static load of over 1 tonne, with an elastic limit of over 250 kg.

One femoral component passed the fatigue test over 10 million cycles under a load of 900 N in identical installation conditions.

### Bibliography:

- [1] Arima J., Whiteside A., Mac Carthy D.S - Femoral rotational alignment based on the anteroposterior axis, in total knee arthroplasty in a valgus knee. J.Bone Joint Surgery Am 1995, 77 : 1331-1334.
- [2] Olcott C.W., Scott R.D. - A comparison of 4 intraoperative methods to determine femoral component rotation during total knee arthroplasty. J. Arthroplasty 2000 ; 15 : 22-26.
- [3] Bellemans J., Robijns F., Duerinckx J., Banks S., Vandenuecker H. The influence of tibial slope on maximal flexion after total knee arthroplasty. Knee Surg. Sports Traumatol. Arthrosc. 2005 ; 13 : 193-196.
- [4] Boisgard S., Moreau Pe., Deschamps S., Courthalc C., Silbert H., Moreel P., Michel J.L., Levai Jp. - Computed tomographic study of the posterior condylar angle in arthritic knees : its use in the rotational positioning of the femoral implant of total knee prosthesis. - Surg. Radiol. Anat. 2003 ; 25, 330-334.
- [5] Dennis D.A., Channer M., Susman H., Stringer E.A. - Intramedullary versus extramedullary tibial alignment systems in total knees arthroplasty. - J. Arthroplasty 1993 ; 8 : 43-47.
- [6] Vertullo C.J., Easley M.E., Scott W.N., Insall J.N. - Mobile bearings in primary knee arthroplasty. - J. Am. Acad. Orthop. Surg. 2001 Nov-Dec ; 9(6) : 355-64.
- [7] Ranawat C.S., Komistek R.D., Rodriguez J.A., Dennis D.A., Anderle M. - In vivo kinematics for fixed and mobile-bearings posterior stabilized knee prosthesis. - Clin. Orthop. 2004 ; 418 : 184-190.
- [8] Morawa L.G., Carpenter C.W., Desai K.S., Kerr D.R., McClure M.G. - Use of the single medial/lateral radius design in cemented total knee arthroplasty : Clinical and radiographic results at 5-year follow-up. - Seminars in Arthroplasty vol. 7 N° 3 (july) 1996 : 276-284.
- [9] Mihalko W.M., Whiteside L.A., Krackow K.A. - Comparison of ligament balancing techniques during total knee arthroplasty. - J. Bone Joint Surg. Am 2003, 85 suppl. 132-135.
- [10] Star M.J., Kaufman K.R., Irbyse Colwell C.W. - The effect of patellar thickness on patella femoral forces after resurfacing. - Clin. Orthop. 1996 ; 322, 279-285.
- [11] Kuster M.S., Horz S., Spalinger E., Stachowiak W., Gächter A. - The effects of conformity and load in total knee replacement. - Cli. Orthop. 2000 ; 375, 302-312
- [12] Li S., Scuderi G., Furman B.D., Bhattacharyya S., Schmiegg J.J., Insall J.N. - Assessment of backside wear from the analysis of 55 retrieved tibial inserts. - Cli. Orthop. 2002 ; 404, 75-82.
- [13] Bartel D.L., Bicknell V.L., Ithaca, Wright T.M. - The effect of conformity, thickness and material on stresses in ultra-high molecular weight components for total joint replacement. - J. Bone and joint surg. 1986.

# SURGICAL TECHNIQUE

The following surgical technique describes the steps of a **MADISON** traditional primary total knee prosthesis procedure.

The instrumentation is designed solely for inserting the **MADISON** total knee prosthesis.

The Manufacturer shall not be liable should a **MADISON** total knee prosthesis procedure be performed without the specific instruments.

## STEPS

### Important:

**a 1.27 mm-thick saw blade is recommended to ensure the quality of bone cuts.**

1. Make a proximal tibial cut of at least 10 mm using the extramedullary or intramedullary alignment device and a right or left cutting guide with a 0 or 4° posterior slope.
2. Make a distal femoral cut of 8 mm (10 or 12 mm if necessary) using the distal cutting guide in left or right position, fixed on the alignment guide previously set to 3 to 9 degrees according to the HKS angle.
3. Measure the external rotation angle with the Whiteside's line angle finder.
4. Choose the AP and ML femoral size and set the external rotation based on the measurement performed in step 3.
5. Drill positioning holes for the 4-in-1 femoral cutting guide.
6. Make femoral cuts using the 4-in-1 guide selected in step 4: anterior cortex, posterior condyles, anterior and posterior chamfers.
7. Cut the PS cage (not for UC procedure) and drill holes for the femoral pegs.

Preoperative planning may be carried out using radiographic templates, enabling the surgeon to evaluate the size of the femoral and tibial components and the length of the tibial extension rod (if this option is considered).

Radiographic templates are available with a 1:1 or 1:1.15 scale.


Incision is performed according to the surgeon's preferred practice.

The operation may start independently with the tibia or the femur. The choice of the tibial intramedullary or extramedullary option also depends on the surgeon's preferred practice.

8. Check the size of the tibial baseplate using the range of test plates (see component compatibility tables if required).
9. Use a CR trial femur for UC procedures, or a PS trial femur for the PS procedure. Drill holes for the femoral pegs if using a CR femur.
10. Conduct trial prosthesis tests after snapping a trial insert onto the corresponding trial plate.
11. Check leg alignment and ligament balance in extension and flexion.
12. Prepare the tibial keel (cemented or cementless).
13. Prepare the patella with an 8 or 10 mm cut.
14. Choose the size and thickness of the patella (8 or 10 mm), drill holes for the three anchor pegs and test the selected patella.
15. Insert the implants (recommended order): tibial baseplate – femoral component – tibial insert – patella.

# INCISION AND EXPOSURE

- A** An internal, vertical cut is made into the skin.



The surgeon has a choice of three types of surgical approach to insert the **MADISON** knee prosthesis:


- B** Traditional internal parapatellar
- C** Midvastus
- D** Subvastus

# TIBIAL PREPARATION

## EXTRAMEDULLARY OPTION

The extramedullary alignment guide is assembled as shown. **Do not forget to install the height adjustment knob on the pole before inserting it into the alignment support.**

Finally, the tibial cutting guide is inserted over the top of the pole.





### Option:

An alignment extension is available as an option for long tibias.

**1** The extramedullary alignment guide is a rigid framework that must be set up according to the following instructions:

- Position the clamps on the two malleoli.
- Position the pole by inserting the long spike into the anterior cruciate ligament insertion (or behind).
- Position the V against the tibial crest.

The positioning of the guide eliminates tibial slope and rotation.



**2** The surgeon chooses a cutting block with a 0° or 4° posterior tibial slope according to their desired final slope.




**Reminder:**

The **MADISON** tibial insert has a 4° built-in slope.

The selected tibial cutting block is inserted into the proximal part of the extramedullary guide.


- 3** The posterior slope is adjusted in relation to the parallelism between the tibia profile axis and the distal alignment rod mounted on the removable arch.

Once the position is defined, the pole is permanently fixed in place by impacting the two spikes into the tibial plateau.



- 4** The cutting height is determined using one of the two styluses on the support and placed on the cutting guide:

- 2 mm stylus on the worn side
- 10 mm stylus on the healthy side




Adjust the stylus height using the height adjustment knob .  
Lock the cutting block in position once the resection level is chosen.  
Fix it into position by inserting two pins into the holes marked with "0".  
It can be repositioned on the cutting guide by +2 or +4 mm.

- 5** After releasing all the knobs on the system and loosening the height adjustment knob, the pole can be extracted using the extractor then manually removed in order to avoid moving the cutting guide.



- 6** The cutting block is then brought into contact with the anterior surface of the tibia and fixed with a third headed pin.




# TIBIAL PREPARATION

## INTRAMEDULLARY OPTION

- 1 The proximal tibia is drilled using the 8 mm bit in the tibial spine in line with the canal. The 8 mm centromedullary rod attached to the T-handle is inserted.

The T-handle is then removed.



2



A cutting guide with a 0° or 4° posterior tibial slope is chosen according to the surgeon's desired slope.

The 8 mm adapter is installed on the support.

The tibial cutting guide is installed on the intramedullary support and the assembly is positioned on the intramedullary rod.

An optional fixing pin inserted on the front part of the support prevents it from rotating in the sagittal plane.


The cutting guide is brought as close to the anterior surface of the tibia as possible.



3

The cutting height is determined using one of the two styluses on the support placed on the cutting guide:



- 2 mm stylus on the worn side
- 10 mm stylus on the healthy side




10 mm stylus

2 mm stylus

- 4** The guide is fixed into position using 2 pins inserted into the holes marked "0".  
If necessary, + 2 mm and + 4 mm holes are available for additional tibial cutting.





- 5** Remove the rod using the T-handle then remove the support.  
The tibial cutting block is brought close to the anterior surface of the tibia then locked using a headed pin.



# TIBIAL PREPARATION

## PROXIMAL CUT

The proximal tibial cut is made following the extramedullary or intramedullary option.





At this stage, the size of the tibial baseplate can be estimated using trial baseplates and by turning over the resected tibial section.

# FEMORAL PREPARATION


## DISTAL CUT

- 1 Set the degree of femoral valgus on the distal alignment guide (HKS angle calculated preoperatively): right or left side 3° to 9° (1° increments).



The 8 or 10 mm distal femoral cut can be selected by placing the cutting guide support in the 8 or 10 mm position on the alignment guide. The recommended cut is 8 mm.

- 2 The distal femoral cutting block is assembled on its support then positioned on the alignment guide.



8 mm cut

10 mm cut


- 3** Drill the entrance hole into the notch using the 8 mm diameter bit, in line with the centromedullary femoral canal.

Insert the 8 mm diameter centromedullary rod using the removable T-handle, ensuring that at least 10 cm is left outside the bone.

Remove the T-handle and insert the distal femoral alignment guide on the rod until it comes into contact with at least one distal femoral condyle.



Fix the cutting block in position by inserting two pins into the "O" holes.

Loosen the knob on the femoral guide support, remove the centromedullary rod using the T-handle then remove the alignment guide.



- 4** The guide is brought as close as possible to the anterior part of the femur then locked with a headed pin.

The distal femoral cut is made then the guide is removed.




# SPACERS OPTION

The Extension (EXT) base spacer represents the minimum space required for a **MADISON** prosthesis, which is 18 mm (8 mm femur + baseplate and 10 mm insert).

The Extension (EXT) base spacer is used to check the space and full extension when the tibial proximal cut and distal femoral cuts are made.

The Flexion (FLEX) base spacer represents the minimum space required for a tibial baseplate and the smallest insert thickness, which is 10 mm.

The Flexion (FLEX) base spacer is used at 90° flexion to check the space between the proximal tibial cut and the unresected posterior condyles.




When all the femoral and tibial cuts have been made, the extension (EXT) base spacer may be used to check the spaces in extension and flexion.

If necessary, a + 2 mm or + 4 mm spacer may be inserted on the flexion and extension base spacers (FLEX and EXT).



# FEMUR PREPARATION

## MEASUREMENT AND EXTERNAL ROTATION

- 1 An angle finder is used to determine the external femoral rotation angle, based on the Whiteside's line principle and in relation to the posterior bicondylar line.



Rotation angle/Whiteside's line.



The calculated angle will be set on the Evo Sizer. In this example, a 3° external rotation angle is read and set for a right knee.


2

The Evo Sizer instrument has 3 components:

- Femoral/rotation Evo Sizer - Base
- Femoral/rotation Evo Sizer - Support
- Femoral/rotation Evo Sizer - Stylus




The stylus tip is flush with the notch.




The stylus tip indicates the saw blade exit point.

1. Position of the Evo Sizer against the surface of the distal femoral cut with the knee in 90° flexion.

2. The Support foot is in contact with the posterior femoral condyles.



The AP size is determined when the stylus comes into contact with the anterior cortex, the line on the stylus support indicates the size in the upper window of the Evo Sizer base.





The ML size is checked by aligning the lateral arms with the bone edges. ML 1, 2, 3, ...8 indicate the ML size of the corresponding implant (ensure that each arm shows the same figure).

### First case


ML = size 4 / AP = size 4

Set the stylus to position 4 and place the stylus + support assembly onto the Evo Sizer Base.





In this case, if size 4 is confirmed with the AP and ML measurement, drill into the central holes using the 3.2 mm bit.

**Next step: 4-in-1 cutting guide, size 4 (page 25).**



### Important:

This only applies to changing from size 3 up to 7 (2 mm increment between sizes).





### Second case

ML = size 4 / AP = size 5

In this case, the AP measurement indicates size 5 but the ML measurement indicates size 4.

If the choice is made to maintain size 4, the holes will be drilled in the upper holes of the Evo Sizer base in order to avoid creating a notch in the anterior cortex.

**Next step: 4-in-1 cutting guide, size 4.**



### Third case


ML = size 4 / AP = size 3

In this case, the AP measurement indicates size 3 but the ML measurement indicates size 4.

If the choice is made to maintain size 4, the holes will be drilled in the lower holes of the Evo Sizer Base in order to prevent the prosthesis being raised in relation to the anterior cortex.

## 4-IN-1 CUTTING GUIDE

- 1 The 4-in-1 cutting guide for the selected size is positioned in the previously created holes.



2

- The guide is pushed into contact with the distal femoral cut.


A final check can be carried out on the anterior cortex using the test plate. Use two headed fixing pins in a convergent position to secure the guide in position.

## ANTERIOR, POSTERIOR AND CHAMFERED CUTS


- 3** Cuts must be made using a saw blade 1.27 mm thick and ideally 25 mm wide and 90 mm long.

**The following order is recommended:**

- A** anterior cut
- B** posterior cuts
- C** anterior chamfer
- D** posterior chamfers



- 4** Remove the headed pins and cutting guide using the universal extractor.




## CR PROCEDURE


### Only for UC

**1** The trial CR femoral implant in the selected size is impacted using the femoral impactor.

**2** The impactor is removed when the ML positioning of the trial implant is determined.






**3** Two holes for the femoral pegs are drilled through the trial implant using the 6 mm drill with stop.



## PS PROCEDURE

### PS cage and implant positioning

- 1 Three headed pins are inserted to lock the PS cage cutting guide into position.
- 2 Start by cutting with a saw blade, using the two internal edges of the guide and not exceeding a depth of 4 mm 1.
- 3 The PS cage chisel corresponding to the size of the guide is inserted into the upper guide slot and gently impacted up to the depth stop to complete PS cage cutting 2.
- 4 The chisel is removed in a straight line using the universal extractor which then enables the three headed pins to be removed.




- 5 Two holes for the femoral pegs are drilled through the PS cage guide using the 6 mm drill with stop.

# TESTS



- 1 The selected trial tibial baseplate is temporarily fixed onto the tibial cut using three headed pins.

Rotation can be checked using the external rod inserted into the handle.



- 2 The trial femoral component is impacted again using the specific impactor.


## Important :

The CR version must be used for the UC procedure.  
The PS version must be used for the PS procedure.

- 3 The trial insert corresponding to the chosen baseplate is snapped into position.

The thinnest trial insert is usually used first (10 mm). Depending on the required degree of stability, a 12 to 20 mm\*\* insert for UC and PS can be used.

The trial insert may be removed using the gripper extractor tool inserted into the 2 slots.




PS


UC

\*\* Some thicknesses are optional. For more details please contact your local representative.

- 4** The centre line of the lower limb and ligament stability in extension and flexion can be checked.



- 5** The definitive position of the trial tibial baseplate is marked on the anterior edge of the tibia using an electric scalpel.



# TIBIA

## KEEL PREPARATION



The insert and trial femur are removed.

The guide pillar is placed on the two trial baseplate pegs.

There is a cemented and cementless version of the tibial keel impactor.

The selected impactor is inserted into the pillar and gently impacted up to the depth stop on the pillar.

Remove the impactor in a straight line using the universal extractor, enabling the three headed pins to be removed.



# EXTENSION ROD OPTION




If a tibial extension rod is required, the specific drilling guide is positioned on the two trial baseplate pegs.

Insert the sleeve with the rod size marking.

Pull the knob out to fully insert the sleeve then slowly lift the sleeve to lock it into the selected rod position.

In this example, we have selected a 35 mm-long rod for a size 3, 4 or 5 tibial baseplate.

Then drill up to the stop using the 11 mm diameter bit.



## Positioning the definitive extension rod:

unscrew the cap under the tibial keel using a 5 mm hex screwdriver and screw in the chosen extension.



# PATELLA

It is the surgeon's choice whether or not to resurface the patella.

The patella clamp is assembled and the resection level is set to 8 or 10 mm.

The clamp is firmly locked before cutting the patella.


The patella must be cut in one go as removing the blade during resection may alter the cutting level (pressure-responsive clamp).



# IMPLANTS

The recommended order for positioning implants is:

- 1 Tibial baseplate
- 2 Femoral component
- 3 Tibial insert
- 4 Patella




Trial and implant femoral impactor

The knee is then placed in full extension which maintains pressure on the the femur and tibial baseplate to ensure correct impaction.

If the patella is in place, cement is applied to the implant or the bone cut according to the surgeon's usual practice.

The implant is then placed in the three holes and pressure is applied using the patella clamp with cement nozzle.

The femoral impactor is used only to insert the trial and definitive implant.



# IMPLANT DIMENSIONS

Unit: mm

## FEMUR

|                          | 1*   | 2    | 3    | 4    | 5    | 6    | 7    | 8    |
|--------------------------|------|------|------|------|------|------|------|------|
| Anteroposterior (AP)     | 53,5 | 57,5 | 61,5 | 63,5 | 65,5 | 67,5 | 70   | 74   |
| Mediolateral (ML)        | 53,5 | 57   | 61,0 | 65,0 | 68,5 | 72,0 | 74,0 | 81,5 |
| Intercondylar notch (IN) | 16,5 | 16,5 | 17,0 | 17,5 | 18,0 | 18,5 | 19,0 | 19,5 |

\* Size 1 is in the Madison range

## TIBIA

|                      | 1    | 2    | 3    | 4    | 5    | 6    | 7    | 8    |
|----------------------|------|------|------|------|------|------|------|------|
| Anteroposterior (AP) | 40,5 | 43,5 | 46,5 | 49,5 | 52,5 | 54,0 | 56,0 | 59,5 |
| Mediolateral (ML)    | 61,5 | 66,0 | 70,5 | 75,0 | 80,0 | 82,5 | 85,5 | 90,5 |

## PS TIBIAL INSERTS

|            | 1  | 2  | 3    | 4  | 5    | 6  | 7    | 8  |
|------------|----|----|------|----|------|----|------|----|
| Thickness  | 10 | 10 | 10   | 10 | 10   | 10 | 10   | 10 |
|            | 12 | 12 | 12   | 12 | 12   | 12 | 12   | 12 |
|            | 14 | 14 | 14   | 14 | 14   | 14 | 14   | 14 |
| PS peg (D) | 15 | 15 | 15,5 | 16 | 16,5 | 17 | 17,5 | 18 |


## UC TIBIAL INSERTS

|               | 1  | 2  | 3  | 4  | 5  | 6  | 7  | 8  |
|---------------|----|----|----|----|----|----|----|----|
| Thickness (T) | 10 | 10 | 10 | 10 | 10 | 10 | 10 | 10 |
|               | 12 | 12 | 12 | 12 | 12 | 12 | 12 | 12 |
|               | 14 | 14 | 14 | 14 | 14 | 14 | 14 | 14 |

\*\* Please check product availability with your local representative.




## PATELLA

| DIAMETER      | 30 | 33 | 36 | 39 |
|---------------|----|----|----|----|
| Peg diameter  | 5  | 5  | 5  | 5  |
| Thickness (T) | 8  | 8  | 8  | 8  |
|               | 10 | 10 | 10 | 10 |



## TIBIAL RODS

| DIAMETER   | 11 |
|------------|----|
| Length (L) | 35 |
|            | 55 |



# IMPLANT SIZES



**Insert size**



## MOBILE CONFIGURATION

FEMORAL COMPONENT SIZE

|   | 1* | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
|---|----|---|---|---|---|---|---|---|
| 1 | 1  | 1 | 2 | 2 |   |   |   |   |
| 2 | 1  | 2 | 2 | 3 | 3 |   |   |   |
| 3 |    | 2 | 3 | 3 | 4 | 4 |   |   |
| 4 |    |   | 3 | 4 | 4 | 5 | 5 |   |
| 5 |    |   |   | 4 | 5 | 5 | 6 | 6 |
| 6 |    |   |   |   | 5 | 6 | 6 | 7 |
| 7 |    |   |   |   |   | 6 | 7 | 7 |
| 8 |    |   |   |   |   |   | 7 | 8 |

\* Size 1 is in the **MADISON** range

### Important:

When you use a tibial baseplate two sizes smaller than the femur, the insert must be one size bigger than the tibial baseplate.



## FIXED CONFIGURATION

### FEMORAL COMPONENT SIZE

|   | 1* | 2 | 3 | 4 | 5 | 6 | 7 | 8 |
|---|----|---|---|---|---|---|---|---|
| 1 | 1  | 1 |   |   |   |   |   |   |
| 2 | 2  | 2 | 2 |   |   |   |   |   |
| 3 |    | 3 | 3 | 3 |   |   |   |   |
| 4 |    |   | 4 | 4 | 4 |   |   |   |
| 5 |    |   |   | 5 | 5 | 5 |   |   |
| 6 |    |   |   |   | 6 | 6 | 6 |   |
| 7 |    |   |   |   |   | 7 | 7 | 7 |
| 8 |    |   |   |   |   |   | 8 | 8 |

\* Size 1 is in the **MADISON** range

The patella components are all compatible with the range of femoral components.




# INSTRUMENTATION



# MADISON EVO

## BOX 1

| NO | QTY | DESCRIPTION                                          |
|----|-----|------------------------------------------------------|
| 1  | 1   | <b>540009</b> - Distal femoral cutting guide         |
| 2  | 1   | <b>540010</b> - Femoral valgisation guide            |
| 3  | 1   | <b>540011</b> - Distal femoral cutting guide support |
| 4  | 1   | <b>540376</b> - 4-in-1 semi-open cutting guide, S2   |
| 5  | 1   | <b>540377</b> - 4-in-1 semi-open cutting guide, S3   |
| 6  | 1   | <b>540378</b> - 4-in-1 semi-open cutting guide, S4   |
| 7  | 1   | <b>540379</b> - 4-in-1 semi-open cutting guide, S5   |
| 8  | 1   | <b>540380</b> - 4-in-1 semi-open cutting guide, S6   |
| 9  | 1   | <b>540381</b> - 4-in-1 semi-open cutting guide, S7   |
| 10 | 1   | <b>540530</b> - Femoral/rotation Evo Sizer - Base    |
| 11 | 1   | <b>540531</b> - Femoral/rotation Evo Sizer - Support |
| 12 | 1   | <b>540532</b> - Femoral/rotation Evo Sizer - Stylus  |
| 13 | 1   | <b>540041</b> - Whiteside's line angle finder        |
| 14 | 1   | <b>540043</b> - Centromedullary rod, dia. 8 mm       |
| 15 | 1   | <b>540044</b> - Hudson femoral drill, 8 mm           |
| 16 | 1   | <b>540396</b> - Hudson femoral drill, 8 mm           |
| 17 | 1   | <b>540398</b> - Hudson femoral peg drill, 6 mm       |




## BOX 2 TOP

### "ML cutting guide" composition

| No | QTY | DESCRIPTION                                     |
|----|-----|-------------------------------------------------|
| 1  | 1   | 540137 - Hex screwdriver, dia. 5 mm             |
| 2  | 1   | 540550 - EVO tibial stylus, 2 mm                |
| 3  | 1   | 540551 - EVO tibial stylus, 10 mm               |
| 4  | 1   | 540112 - Tibial intramedullary alignment        |
| 5  | 1   | 540113 - Height adjustment knob                 |
| 6  | 1   | 540115 - Rod adaptor, dia. 8 mm                 |
| 7  | 1   | 540124 - Extramedullary pole                    |
| 8  | 1   | 540125 - Extramedullary alignment support       |
| 9  | 1   | 540126 - V-shaped extramedullary positioner     |
| 10 | 1   | 540127 - Short distal tibial alignment rod      |
| 11 | 1   | 540128 - Removable distal tibial alignment arch |
| 12 | 1   | 540129 - V-shaped tibial positioner stop        |
| 13 | 1   | 540130 - Malleolar clamp, part 1/2              |
| 14 | 1   | 540131 - Malleolar clamp, part 2/2              |
| 15 | 1   | 540132 - Malleolar clamp stop                   |
| 16 | 1   | 540173 - Fixing pin extractor                   |
| 17 | 1   | 540403 - ML cutting guide, right, 0°            |
| 18 | 1   | 540404 - ML cutting guide, left, 0°             |
| 19 | 1   | 540405 - ML cutting guide, right, 4°            |
| 20 | 1   | 540406 - ML cutting guide, left, 4°             |
| 21 | 1   | 540549 - EVO tibial stylus support              |
| 22 | 1   | 540392 - Quick pin connector                    |


### "HEMI GUIDE" composition

| No | QTY | DESCRIPTION                                     |
|----|-----|-------------------------------------------------|
| 1  | 1   | 540137 - Hex screwdriver, dia. 5 mm             |
| 2  | 1   | 540104 - Reference tibial stylus, 2 mm          |
| 3  | 1   | 540105 - Reference tibial stylus, 10 mm         |
| 4  | 1   | 540112 - Tibial intramedullary alignment        |
| 5  | 1   | 540113 - Height adjustment knob                 |
| 6  | 1   | 540115 - Rod adaptor, dia. 8 mm                 |
| 7  | 1   | 540124 - Extramedullary pole                    |
| 8  | 1   | 540125 - Extramedullary alignment support       |
| 9  | 1   | 540126 - V-shaped extramedullary positioner     |
| 10 | 1   | 540127 - Short distal tibial alignment rod      |
| 11 | 1   | 540128 - Removable distal tibial alignment arch |
| 12 | 1   | 540129 - V-shaped tibial positioner stop        |
| 13 | 1   | 540130 - Malleolar clamp, part 1/2              |
| 14 | 1   | 540131 - Malleolar clamp, part 2/2              |
| 15 | 1   | 540132 - Malleolar clamp stop                   |
| 16 | 1   | 540173 - Fixing pin extractor                   |
| 17 | 1   | 540106 - Tibial cutting guide, left, 4°         |
| 18 | 1   | 540108 - Tibial cutting guide, right, 0°        |
| 19 | 1   | 540109 - Tibial cutting guide, left, 0°         |
| 20 | 1   | 540110 - Tibial cutting guide, right, 4°        |
| 21 | 1   | 540111 - Tibial stylus guide                    |
| 22 | 1   | 540392 - Quick pin connect                      |




## BOX 2 BOTTOM

| No | QTY | DESCRIPTION                               |
|----|-----|-------------------------------------------|
| 1  | 1   | 540002 - Trial tibial baseplate, S2       |
| 2  | 1   | 540003 - Trial tibial baseplate, S3       |
| 3  | 1   | 540004 - Trial tibial baseplate, S4       |
| 4  | 1   | 540005 - Trial tibial baseplate, S5       |
| 5  | 1   | 540006 - Trial tibial baseplate, S6       |
| 6  | 1   | 540007 - Trial tibial baseplate, S7       |
| 7  | 1   | 540107 - Resection test plate             |
| 8  | 1   | 540116 - Tibial keel impactor pillar      |
| 9  | 1   | 540117 - Press-fit pin impactor, S1-2     |
| 10 | 1   | 540118 - Press-fit pin impactor, S3-4-5   |
| 11 | 1   | 540119 - Press-fit pin impactor, S6-7-8   |
| 12 | 1   | 540120 - Cemented pin impactor, S1-2      |
| 13 | 1   | 540121 - Cemented pin impactor, S3-4-5    |
| 14 | 1   | 540122 - Cemented pin impactor, S6-7-8    |
| 15 | 1   | 540370 - Tibial trial gripper/extractor   |
| 16 | 1   | 540138 - Tibial drill bushing             |
| 17 | 1   | 540139 - Tibial drill bushing base        |
| 18 | 1   | 540397 - Hudson femoral drill, 11 mm      |
| 19 | 1   | 540141 - External alignment rod 1         |
| 20 | 1   | 540282 - External alignment rod 2         |
| 21 | 1   | 540283 - Pin box                          |
|    | 1   | 540284 - Pin box cover                    |
|    | 3   | 540162 - Fixing nail, dia. 3.2 mm, L1     |
|    | 3   | 540163 - Fixing nail, dia. 3.2 mm, L2     |
|    | 3   | 540165 - Fixing pin, dia. 3.2 mm, L2      |
|    | 3   | 540393 - Quick connect twist pin, L1      |
|    | 3   | 540394 - Quick connect twist pin, L2      |
|    | 3   | 540395 - Quick connect twist pin, L3      |
|    | 1   | 540402 - Drill bit, 3.2 mm, AO connection |
|    |     |                                           |




## BOX 3 TOP

| No | QTY | DESCRIPTION                          |
|----|-----|--------------------------------------|
| 1  | 1   | 540036 - PS cage chisel, N°1         |
| 2  | 1   | 540037 - PS cage chisel, N°2         |
| 3  | 1   | 540038 - PS cage chisel, N°3         |
| 4  | 1   | 540211 - PS trial insert, S2 - 10 mm |
| 5  | 1   | 540212 - PS trial insert, S2 - 12 mm |
| 6  | 1   | 540213 - PS trial insert, S2 - 14 mm |
| 7  | 1   | 540214 - PS trial insert, S3 - 10 mm |
| 8  | 1   | 540215 - PS trial insert, S3 - 12 mm |
| 9  | 1   | 540216 - PS trial insert, S3 - 14 mm |
| 10 | 1   | 540217 - PS trial insert, S4 - 10 mm |
| 11 | 1   | 540218 - PS trial insert, S4 - 12 mm |
| 12 | 1   | 540219 - PS trial insert, S4 - 14 mm |
| 13 | 1   | 540220 - PS trial insert, S5 - 10 mm |
| 14 | 1   | 540221 - PS trial insert, S5 - 12 mm |
| 15 | 1   | 540222 - PS trial insert, S5 - 14 mm |
| 16 | 1   | 540223 - PS trial insert, S6 - 10 mm |
| 17 | 1   | 540224 - PS trial insert, S6 - 12 mm |
| 18 | 1   | 540225 - PS trial insert, S6 - 14 mm |
| 19 | 1   | 540226 - PS trial insert, S7 - 10 mm |
| 20 | 1   | 540227 - PS trial insert, S7 - 12 mm |
| 21 | 1   | 540228 - PS trial insert, S7 - 14 mm |




## BOX 3 BOTTOM

| No | QTY | DESCRIPTION                            |
|----|-----|----------------------------------------|
| 1  | 1   | 540029 - PS cage cutting guide, S2     |
| 2  | 1   | 540030 - PS cage cutting guide, S3     |
| 3  | 1   | 540031 - PS cage cutting guide, S4     |
| 4  | 1   | 540032 - PS cage cutting guide, S5     |
| 5  | 1   | 540033 - PS cage cutting guide, S6     |
| 6  | 1   | 540034 - PS cage cutting guide, S7     |
| 7  | 1   | 540500 - PS trial EVO femur, S2, Right |
| 8  | 1   | 540501 - PS trial EVO femur, S3, Right |
| 9  | 1   | 540502 - PS trial EVO femur, S4, Right |
| 10 | 1   | 540503 - PS trial EVO femur, S5, Right |
| 11 | 1   | 540504 - PS trial EVO femur, S6, Right |
| 12 | 1   | 540505 - PS trial EVO femur, S7, Right |
| 13 | 1   | 540507 - PS trial EVO femur, S2, Left  |
| 14 | 1   | 540508 - PS trial EVO femur, S3, Left  |
| 15 | 1   | 540509 - PS trial EVO femur, S4, Left  |
| 16 | 1   | 540510 - PS trial EVO femur, S5, Left  |
| 17 | 1   | 540511 - PS trial EVO femur, S6, Left  |
| 18 | 1   | 540512 - PS trial EVO femur, S7, Left  |




## BOX 4 TOP

| No | QTY | DESCRIPTION                           |
|----|-----|---------------------------------------|
| 1  | 1   | 540146 - Trial/drilling end fitting 1 |
| 2  | 1   | 540147 - Trial/drilling end fitting 2 |
| 3  | 1   | 540148 - Trial/drilling end fitting 3 |
| 4  | 1   | 540149 - Trial/drilling end fitting 4 |
| 5  | 1   | 540256 - Trial patella, 30/8 mm       |
| 6  | 1   | 540257 - Trial patella, 30/10 mm      |
| 7  | 1   | 540258 - Trial patella, 33/8 mm       |
| 8  | 1   | 540259 - Trial patella, 33/10 mm      |
| 9  | 1   | 540260 - Trial patella, 36/8 mm       |
| 10 | 1   | 540261 - Trial patella, 36/10 mm      |
| 11 | 1   | 540262 - Trial patella, 39/8 mm       |
| 12 | 1   | 540263 - Trial patella, 39/10 mm      |
| 13 | 1   | 540361 - Flex-ext base plate          |
| 14 | 1   | 540362 - Spacer, 2 mm                 |
| 15 | 1   | 540363 - Spacer, 4 mm                 |
| 16 | 1   | 540371 - Alignment rod sleeve         |
| 17 | 1   | 540418 - Hudson adaptor mandrel       |
| 18 | 1   | 540960 - Townley calliper             |




## BOX 4 BOTTOM

| No | QTY | DESCRIPTION                                   |
|----|-----|-----------------------------------------------|
| 1  | 1   | 540172 - Universal extractor                  |
| 2  | 1   | 540142 - Patella clamp                        |
| 3  | 1   | 540143 - Plate with central pin               |
| 4  | 1   | 540150 - Plate for patella cementing          |
| 5  | 1   | 540153 - Patella plate with pins              |
| 6  | 1   | 540399 - Hudson patella peg drill, 5 mm       |
| 7  | 1   | 540158 - Resection block, 8 mm                |
| 8  | 1   | 540159 - Resection block, 10 mm               |
| 9  | 1   | 540135 - Femoral impactor                     |
| 10 | 1   | 540133 - Baseplate impactor and tibial insert |
| 11 | 1   | 540045 - Trial and implant femoral impactor   |
| 12 | 1   | 540134 - Universal impactor                   |




## BOX 5

| No | QTY | DESCRIPTION                            |
|----|-----|----------------------------------------|
| 1  | 1   | 540514 - CR trial EVO femur, S2, Right |
| 2  | 1   | 540515 - CR trial EVO femur, S3, Right |
| 3  | 1   | 540516 - CR trial EVO femur, S4, Right |
| 4  | 1   | 540517 - CR trial EVO femur, S5, Right |
| 5  | 1   | 540518 - CR trial EVO femur, S6, Right |
| 6  | 1   | 540519 - CR trial EVO femur, S7, Right |
| 7  | 1   | 540521 - CR trial EVO femur, S2, Left  |
| 8  | 1   | 540522 - CR trial EVO femur, S3, Left  |
| 9  | 1   | 540523 - CR trial EVO femur, S4, Left  |
| 10 | 1   | 540524 - CR trial EVO femur, S5, Left  |
| 11 | 1   | 540525 - CR trial EVO femur, S6, Left  |
| 12 | 1   | 540526 - CR trial EVO femur, S7, Left  |
| 13 | 1   | 540288 - UC trial insert, S2 - 10 mm   |
| 14 | 1   | 540289 - UC trial insert, S2 - 12 mm   |
| 15 | 1   | 540290 - UC trial insert, S2 - 14 mm   |
| 16 | 1   | 540291 - UC trial insert, S3 - 10 mm   |
| 17 | 1   | 540292 - UC trial insert, S3 - 12 mm   |
| 18 | 1   | 540293 - UC trial insert, S3 - 14 mm   |
| 19 | 1   | 540294 - UC trial insert, S4 - 10 mm   |
| 20 | 1   | 540295 - UC trial insert, S4 - 12 mm   |
| 21 | 1   | 540296 - UC trial insert, S4 - 14 mm   |
| 22 | 1   | 540297 - UC trial insert, S5 - 10 mm   |
| 23 | 1   | 540298 - UC trial insert, S5 - 12 mm   |
| 24 | 1   | 540299 - UC trial insert, S5 - 14 mm   |
| 25 | 1   | 540300 - UC trial insert, S6 - 10 mm   |
| 26 | 1   | 540301 - UC trial insert, S6 - 12 mm   |
| 27 | 1   | 540302 - UC trial insert, S6 - 14 mm   |
| 28 | 1   | 540303 - UC trial insert, S7 - 10 mm   |
| 29 | 1   | 540304 - UC trial insert, S7 - 12 mm   |
| 30 | 1   | 540305 - UC trial insert, S7 - 14 mm   |




## BOX 6 TOP

| No | QTY | DESCRIPTION                            |
|----|-----|----------------------------------------|
| 1  | 1   | 540208 - PS trial insert, S1, t. 10 mm |
| 2  | 1   | 540209 - PS trial insert, S1, t. 12 mm |
| 3  | 1   | 540210 - PS trial insert, S1, t. 14 mm |
| 4  | 1   | 540229 - PS trial insert, S8, t. 10 mm |
| 5  | 1   | 540230 - PS trial insert, S8, t. 12 mm |
| 6  | 1   | 540231 - PS trial insert, S8, t. 14 mm |
| 7  | 1   | 540285 - UC trial insert, S1, t. 10 mm |
| 8  | 1   | 540286 - UC trial insert, S1, t. 12 mm |
| 9  | 1   | 540287 - UC trial insert, S1, t. 14 mm |
| 10 | 1   | 540306 - UC trial insert, S8, t. 10 mm |
| 11 | 1   | 540307 - UC trial insert, S8, t. 12 mm |
| 12 | 1   | 540308 - UC trial insert, S8, t. 14 mm |
| 13 | 1   | 540001 - Trial tibial baseplate, S1    |
| 14 | 1   | 540008 - Trial tibial baseplate, S8    |



## BOX 6 BOTTOM

| No | QTY | DESCRIPTION                                     |
|----|-----|-------------------------------------------------|
| 1  | 1   | 540375 - EL semi-open femoral cutting guide, S1 |
| 2  | 1   | 540184 - PS trial femur, S1, Left               |
| 3  | 1   | 540176 - PS trial femur, S1, Right              |
| 4  | 1   | 540513 - PS trial EVO femur, S8, Left           |
| 5  | 1   | 540506 - PS trial EVO femur, S8, Right          |
| 6  | 1   | 540200 - CR trial femur, S1, Left               |
| 7  | 1   | 540192 - CR trial femur, S1, Right              |
| 8  | 1   | 540527 - CR trial EVO femur, S8, Left           |
| 9  | 1   | 540520 - CR trial EVO femur, S8, Right          |
| 10 | 1   | 540028 - Cage cutting block, S1                 |
| 11 | 1   | 540035 - Cage cutting block, S8                 |
| 12 | 1   | 540382 - EL semi-open femoral cutting guide, S8 |



# ADDITIONAL EQUIPMENT

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

**540610** Fixed insert impactor

## Box 1 - Femur preparation

**540168** IM drill - dia. 8mm

**540171** 6mm femoral peg drill w.stop

**540600** Curved osteotome

## Box 2 - Tibia preparation

**540140** Tibial preparation drill

**540175** Pin and nail motor connector

**540367** Twist pin, L1

**540368** Twist pain, L2

**540434** Varus-valgus tibial alignment

**540435** IM tibial alignment slope

**540436** Quick tibial knob

**540407** Quick tibial knob

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

## Box 3 - PS trial

**540312** PS trial insert, S2 – 16 mm

**540313** PS trial insert, S2 – 18 mm

**540314** PS trial insert, S2 – 20 mm

**540315** PS trial insert, S3 – 16 mm

**540316** PS trial insert, S3 – 18 mm

**540317** PS trial insert, S3 – 20 mm

**540318** PS trial insert, S4 – 16 mm

**540319** PS trial insert, S4 – 18 mm

**540320** PS trial insert, S4 – 20 mm

**540321** PS trial insert, S5 – 16 mm

**540322** PS trial insert, S5 – 18 mm

**540323** PS trial insert, S5 – 20 mm

**540324** PS trial insert, S6 – 16 mm

**540325** PS trial insert, S6 – 18 mm

**540326** PS trial insert, S6 – 20 mm

**540327** PS trial insert, S7 – 16 mm

**540328** PS trial insert, S7 – 18 mm

**540329** PS trial insert, S7 – 20 mm

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

#### Box 4 - Patella

|               |                                        |
|---------------|----------------------------------------|
| <b>540156</b> | Patella peg drill with stop, dia. 5 mm |
| <b>540364</b> | Spacer, 6 mm                           |
| <b>540365</b> | Spacer, 8 mm                           |
| <b>540366</b> | Spacer, 10 mm                          |
| <b>540409</b> | Tibial extractor                       |
| <b>540452</b> | Rotating patella end fitting           |
| <b>540601</b> | 4-in-1 cutting guide impactor          |

#### Box 5 - UC trials

|               |                             |
|---------------|-----------------------------|
| <b>540336</b> | UC trial insert, S2 - 16 mm |
| <b>540337</b> | UC trial insert, S2 - 18 mm |
| <b>540338</b> | UC trial insert, S2 - 20 mm |
| <b>540339</b> | UC trial insert, S3 - 16 mm |
| <b>540340</b> | UC trial insert, S3 - 18 mm |
| <b>540341</b> | UC trial insert, S3 - 20 mm |
| <b>540342</b> | UC trial insert, S4 - 16 mm |
| <b>540343</b> | UC trial insert, S4 - 18 mm |
| <b>540344</b> | UC trial insert, S4 - 20 mm |
| <b>540345</b> | UC trial insert, S5 - 16 mm |
| <b>540346</b> | UC trial insert, S5 - 18 mm |
| <b>540347</b> | UC trial insert, S5 - 20 mm |
| <b>540348</b> | UC trial insert, S6 - 16 mm |
| <b>540349</b> | UC trial insert, S6 - 18 mm |
| <b>540350</b> | UC trial insert, S6 - 20 mm |
| <b>540351</b> | UC trial insert, S7 - 16 mm |
| <b>540352</b> | UC trial insert, S7 - 18 mm |
| <b>540353</b> | UC trial insert, S7 - 20 mm |

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

#### Box 6 - PS and UC trials 1 and 8

|               |                               |
|---------------|-------------------------------|
| <b>540309</b> | PS trial insert, S1, t. 16 mm |
| <b>540310</b> | PS trial insert, S1, t. 18 mm |
| <b>540311</b> | PS trial insert, S1, t. 20 mm |
| <b>540330</b> | PS trial insert, S8, t. 16 mm |
| <b>540331</b> | PS trial insert, S8, t. 18 mm |
| <b>540332</b> | PS trial insert, S8, t. 20 mm |
| <b>540333</b> | UC trial insert, S1, t. 16 mm |
| <b>540334</b> | UC trial insert, S1, t. 18 mm |
| <b>540335</b> | UC trial insert, S1, t. 20 mm |
| <b>540354</b> | UC trial insert, S8, t. 16 mm |
| <b>540355</b> | UC trial insert, S8, t. 18 mm |
| <b>540356</b> | UC trial insert, S8, t. 20 mm |

# IMPLANTS



## FEMORAL COMPONENTS

### REF. DESCRIPTION

#### Cemented

**140000** Cemented PS femur, S1 L

**140300** Cemented PS EVO femur, S 2 L

**140301** Cemented PS EVO femur, S 3 L

**140302** Cemented PS EVO femur, S 4 L

**140303** Cemented PS EVO femur, S 5 L

**140304** Cemented PS EVO femur, S 6 L

**140305** Cemented PS EVO femur, S 7 L

**140306** Cemented PS EVO femur, S 8 L

#### Cementless

**140016** HA-coated cementless PS femur, S1 L

**140314** HA-coated cementless PS EVO femur, S 2 L

**140315** HA-coated cementless PS EVO femur, S 3 L

**140316** HA-coated cementless PS EVO femur, S 4 L

**140317** HA-coated cementless PS EVO femur, S 5 L

**140318** HA-coated cementless PS EVO femur, S 6 L

**140319** HA-coated cementless PS EVO femur, S 7 L

**140320** HA-coated cementless PS EVO femur, S 8 L

LEFT

### REF. DESCRIPTION

#### Cemented

**140008** Cemented PS femur, S1 R

**140307** Cemented PS EVO femur, S 2 R

**140308** Cemented PS EVO femur, S 3 R

**140309** Cemented PS EVO femur, S 4 R

**140310** Cemented PS EVO femur, S 5 R

**140311** Cemented PS EVO femur, S 6 R

**140312** Cemented PS EVO femur, S 7 R

**140313** Cemented PS EVO femur, S 8 R

#### Cementless

**140024** HA-coated cementless PS femur, S1 R

**140321** HA-coated cementless PS EVO femur, S 2 R

**140322** HA-coated cementless PS EVO femur, S 3 R

**140323** HA-coated cementless PS EVO femur, S 4 R

**140324** HA-coated cementless PS EVO femur, S 5 R

**140325** HA-coated cementless PS EVO femur, S 6 R

**140326** HA-coated cementless PS EVO femur, S 7 R

**140327** HA-coated cementless PS EVO femur, S 8 R

RIGHT



## FEMORAL COMPONENTS

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

### Cemented

**140048** Cemented CR femur, S1 L

**140328** Cemented CR EVO femur, S 2 L

**140329** Cemented CR EVO femur, S 3 L

**140330** Cemented CR EVO femur, S 4 L

**140331** Cemented CR EVO femur, S 5 L

**140332** Cemented CR EVO femur, S 6 L

**140333** Cemented CR EVO femur, S 7 L

**140334** Cemented CR EVO femur, S 8 L

### Cementless

**140064** HA-coated cementless CR femur, S1 L

**140342** HA-coated cementless CR EVO femur, S 2 L

**140343** HA-coated cementless CR EVO femur, S 3 L

**140344** HA-coated cementless CR EVO femur, S 4 L

**140345** HA-coated cementless CR EVO femur, S 5 L

**140346** HA-coated cementless CR EVO femur, S 6 L

**140347** HA-coated cementless CR EVO femur, S 7 L

**140348** HA-coated cementless CR EVO femur, S 8 L

LEFT

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

### Cemented

**140056** Cemented CR femur, S1 R

**140335** Cemented CR EVO femur, S 2 R

**140336** Cemented CR EVO femur, S 3 R

**140337** Cemented CR EVO femur, S 4 R

**140338** Cemented CR EVO femur, S 5 R

**140339** Cemented CR EVO femur, S 6 R

**140340** Cemented CR EVO femur, S 7 R

**140341** Cemented CR EVO femur, S 8 R

### Cementless

**140072** HA-coated cementless CR femur, S1 R

**140349** HA-coated cementless CR EVO femur, S 2 R

**140350** HA-coated cementless CR EVO femur, S 3 R

**140351** HA-coated cementless CR EVO femur, S 4 R


**140352** HA-coated cementless CR EVO femur, S 5 R

**140353** HA-coated cementless CR EVO femur, S 6 R

**140354** HA-coated cementless CR EVO femur, S 7 R

**140355** HA-coated cementless CR EVO femur, S 8 R

RIGHT



## FIXED-BEARING TIBIAL BASEPLATE

| REF.              | DESCRIPTION                                             |
|-------------------|---------------------------------------------------------|
| <b>Cemented</b>   |                                                         |
| 140096            | Cemented fixed-bearing tibial baseplate, S1             |
| 140097            | Cemented fixed-bearing tibial baseplate, S2             |
| 140098            | Cemented fixed-bearing tibial baseplate, S3             |
| 140099            | Cemented fixed-bearing tibial baseplate, S4             |
| 140100            | Cemented fixed-bearing tibial baseplate, S5             |
| 140101            | Cemented fixed-bearing tibial baseplate, S6             |
| 140102            | Cemented fixed-bearing tibial baseplate, S7             |
| 140103            | Cemented fixed-bearing tibial baseplate, S8             |
| <b>Cementless</b> |                                                         |
| 140112            | HA-coated cementless fixed-bearing tibial baseplate, S1 |
| 140113            | HA-coated cementless fixed-bearing tibial baseplate, S2 |
| 140114            | HA-coated cementless fixed-bearing tibial baseplate, S3 |
| 140115            | HA-coated cementless fixed-bearing tibial baseplate, S4 |
| 140116            | HA-coated cementless fixed-bearing tibial baseplate, S5 |
| 140117            | HA-coated cementless fixed-bearing tibial baseplate, S6 |
| 140118            | Embase tib. fixe ss ciment HA T7                        |
| 140119            | Embase tib. fixe ss ciment HA T8                        |



## MOBILE-BEARING TIBIAL BASEPLATE

| REF.              | DESCRIPTION                                              |
|-------------------|----------------------------------------------------------|
| <b>Cemented</b>   |                                                          |
| 140104            | Cemented mobile-bearing tibial baseplate, S1             |
| 140105            | Cemented mobile-bearing tibial baseplate, S2             |
| 140106            | Cemented mobile-bearing tibial baseplate, S3             |
| 140107            | Cemented mobile-bearing tibial baseplate, S4             |
| 140108            | Cemented mobile-bearing tibial baseplate, S5             |
| 140109            | Cemented mobile-bearing tibial baseplate, S6             |
| 140110            | Cemented mobile-bearing tibial baseplate, S7             |
| 140111            | Cemented mobile-bearing tibial baseplate, S8             |
| <b>Cementless</b> |                                                          |
| 140120            | HA-coated cementless mobile-bearing tibial baseplate, S1 |
| 140121            | HA-coated cementless mobile-bearing tibial baseplate, S2 |
| 140122            | HA-coated cementless fixed-bearing tibial baseplate, S3  |
| 140123            | HA-coated cementless mobile-bearing tibial baseplate, S4 |
| 140124            | HA-coated cementless mobile-bearing tibial baseplate, S5 |
| 140125            | HA-coated cementless fixed-bearing tibial baseplate, S6  |
| 140126            | HA-coated cementless mobile-bearing tibial baseplate, S7 |
| 140127            | HA-coated cementless mobile-bearing tibial baseplate, S8 |



## PS INSERTS

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

### Thickness 10 mm

**140155** Universal PS insert, S 1, t. 10 mm

**140158** Universal PS insert, S 2, t. 10 mm

**140161** Universal PS insert, S 3, t. 10 mm

**140164** Universal PS insert, S 4, t. 10 mm

**140167** Universal PS insert, S 5, t. 10 mm

**140170** Universal PS insert, S 6, t. 10 mm

**140173** Universal PS insert, S 7, t. 10 mm

**140176** Universal PS insert, S 8, t. 10 mm

### Thickness 12 mm

**140156** Universal PS insert, S 1, t. 12 mm

**140159** Universal PS insert, S 2, t. 12 mm

**140162** Universal PS insert, S 3, t. 12 mm

**140165** Universal PS insert, S 4, t. 12 mm

**140168** Universal PS insert, S 5, t. 12 mm

**140171** Universal PS insert, S 6, t. 12 mm

**140174** Universal PS insert, S 7, t. 12 mm

**140177** Universal PS insert, S 8, t. 12 mm

### Thickness 14 mm

**140157** Universal PS insert, S 1, t. 14 mm

**140160** Universal PS insert, S 2, t. 14 mm

**140163** Universal PS insert, S 3, t. 14 mm

**140166** Universal PS insert, S 4, t. 14 mm

**140169** Universal PS insert, S 5, t. 14 mm

**140172** Universal PS insert, S 6, t. 14 mm

**140175** Universal PS insert, S 7, t. 14 mm

**140178** Universal PS insert, S 8, t. 14 mm



## UC INSERTS

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

### Thickness 10 mm

**140203** Universal UC insert, S 1, t. 10 mm

**140206** Universal UC insert, S 2, t. 10 mm

**140209** Universal UC insert, S 3, t. 10 mm

**140212** Universal UC insert, S 4, t. 10 mm

**140215** Universal UC insert, S 5, t. 10 mm

**140218** Universal UC insert, S 6, t. 10 mm

**140221** Universal UC insert, S 7, t. 10 mm

**140224** Universal UC insert, S 8, t. 10 mm

### Thickness 12 mm

**140204** Universal UC insert, S 1, t. 12 mm

**140207** Universal UC insert, S 2, t. 12 mm

**140210** Universal UC insert, S 3, t. 12 mm

**140213** Universal UC insert, S 4, t. 12 mm

**140216** Universal UC insert, S 5, t. 12 mm

**140219** Universal UC insert, S 6, t. 12 mm

**140222** Universal UC insert, S 7, t. 12 mm

**140225** Universal UC insert, S 8, t. 12 mm

### Thickness 14 mm

**140205** Universal UC insert, S 1, t. 14 mm

**140208** Universal UC insert, S 2, t. 14 mm

**140211** Universal UC insert, S 3, t. 14 mm

**140214** Universal UC insert, S 4, t. 14 mm

**140217** Universal UC insert, S 5, t. 14 mm

**140220** Universal UC insert, S 6, t. 14 mm

**140223** Universal UC insert, S 7, t. 14 mm

**140226** Universal UC insert, S 8, t. 14 mm



## PATELLA COMPONENTS

| REF.   | DESCRIPTION              |
|--------|--------------------------|
| 140144 | Cemented patella, Ø30/10 |
| 140145 | Cemented patella, Ø33/10 |
| 140146 | Cemented patella, Ø36/10 |
| 140147 | Cemented patella, Ø39/10 |
| 140148 | Cemented patella, Ø30/8  |
| 140149 | Cemented patella, Ø33/8  |
| 140150 | Cemented patella, Ø36/8  |
| 140151 | Cemented patella, Ø39/8  |



## TIBIAL EXTENSION RODS\*

| REF.   | DESCRIPTION               |
|--------|---------------------------|
| 140153 | Tibial rod, Ø11 mm, l. 35 |
| 140154 | Tibial rod, Ø11 mm, l. 55 |



## TIBIAL EXTENSION RODS\* / option

| REF.   | DESCRIPTION                  |
|--------|------------------------------|
| 140276 | Tibial rod, Ø9 mm, l. 35 mm  |
| 140277 | Tibial rod, Ø9 mm, l. 55 mm  |
| 140278 | Tibial rod, Ø9 mm, l. 95 mm  |
| 140275 | Tibial rod, Ø11 mm, l. 95 mm |
| 140279 | Tibial rod, Ø13 mm, l. 35 mm |
| 140280 | Tibial rod, Ø13 mm, l. 55 mm |
| 140281 | Tibial rod, Ø13 mm, l. 95 mm |



## PS INSERTS\*\* / option

| REF.                   | DESCRIPTION                        |
|------------------------|------------------------------------|
| <b>Thickness 16 mm</b> |                                    |
| 140227                 | Universal PS insert, S 1, t. 16 mm |
| 140230                 | Universal PS insert, S 2, t. 16 mm |
| 140233                 | Universal PS insert, S 3, t. 16 mm |
| 140236                 | Universal PS insert, S 4, t. 16 mm |
| 140239                 | Universal PS insert, S 5, t. 16 mm |
| 140242                 | Universal PS insert, S 6, t. 16 mm |
| 140245                 | Universal PS insert, S 7, t. 16 mm |
| <b>Thickness 18 mm</b> |                                    |
| 140228                 | Universal PS insert, S 1, t. 18 mm |
| 140231                 | Universal PS insert, S 2, t. 18 mm |
| 140234                 | Universal PS insert, S 3, t. 18 mm |
| 140237                 | Universal PS insert, S 4, t. 18 mm |
| 140240                 | Universal PS insert, S 5, t. 18 mm |
| 140243                 | Universal PS insert, S 6, t. 18 mm |
| 140246                 | Universal PS insert, S 7, t. 18 mm |
| <b>Thickness 20 mm</b> |                                    |
| 140229                 | Universal PS insert, S 1, t. 20 mm |
| 140232                 | Universal PS insert, S 2, t. 20 mm |
| 140235                 | Universal PS insert, S 3, t. 20 mm |
| 140238                 | Universal PS insert, S 4, t. 20 mm |
| 140241                 | Universal PS insert, S 5, t. 20 mm |
| 140244                 | Universal PS insert, S 6, t. 20 mm |
| 140247                 | Universal PS insert, S 7, t. 20 mm |

\* Please check product availability with your local representative.

\*\* Some thicknesses are optional. For more details please contact your local representative.



## UC\*\* INSERTS / option

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

### Thickness 16 mm

|               |                                    |
|---------------|------------------------------------|
| <b>140251</b> | Universal UC insert, S 1, t. 16 mm |
| <b>140254</b> | Universal UC insert, S 2, t. 16 mm |
| <b>140257</b> | Universal UC insert, S 3, t. 16 mm |
| <b>140260</b> | Universal UC insert, S 4, t. 16 mm |
| <b>140263</b> | Universal UC insert, S 5, t. 16 mm |
| <b>140266</b> | Universal UC insert, S 6, t. 16 mm |
| <b>140269</b> | Universal UC insert, S 7, t. 16 mm |

### Thickness 18 mm

|               |                                    |
|---------------|------------------------------------|
| <b>140252</b> | Universal UC insert, S 1, t. 18 mm |
| <b>140255</b> | Universal UC insert, S 2, t. 18 mm |
| <b>140258</b> | Universal UC insert, S 3, t. 18 mm |
| <b>140261</b> | Universal UC insert, S 4, t. 18 mm |
| <b>140264</b> | Universal UC insert, S 5, t. 18 mm |
| <b>140267</b> | Universal UC insert, S 6, t. 18 mm |
| <b>140270</b> | Universal UC insert, S 7, t. 18 mm |

| REF. | DESCRIPTION |
|------|-------------|
|------|-------------|

### Thickness 20 mm

|               |                                    |
|---------------|------------------------------------|
| <b>140253</b> | Universal UC insert, S 1, t. 20 mm |
| <b>140256</b> | Universal UC insert, S 2, t. 20 mm |
| <b>140259</b> | Universal UC insert, S 3, t. 20 mm |
| <b>140262</b> | Universal UC insert, S 4, t. 20 mm |
| <b>140265</b> | Universal UC insert, S 5, t. 20 mm |
| <b>140268</b> | Universal UC insert, S 6, t. 20 mm |
| <b>140271</b> | Universal UC insert, S 7, t. 20 mm |

\*\* Some thicknesses are optional. For more details please contact your local representative.

# NOTES

[illegible]

# NOTES

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

All the medical devices mentioned in this document are CE marked according to Directive 93/42/EEC and its amendments, unless specifically identified as "not CE marked".

The medical devices mentioned in this document are class III, IIb, IIa, and I devices. Class III, IIb, and IIa medical devices are marked CE 0459 by GMED.

Before using any product, please read the instructions for use and surgical technique. Refer to the labels and instruction leaflets for the complete list of indications, contraindications, risks, warnings, precautions and instructions for use. Reimbursement status in France: the implants are included on the List of Reimbursable Products and Services (Liste des Produits et Prestations Remboursables, LPPR) referred to in article L165-1 of the Social Security Code and are thus eligible for reimbursement by the health insurance scheme.

©2021 SERF. All rights reserved.

DISTRIBUTED BY

CE 0459



IMPLANET  
Technopole Bordeaux Montesquieu  
Allée François Magendie • 33650 Martillac - FRANCE

DISTRIBUTED BY SERF  
85 avenue des Bruyères, 69150 Décines-Charpieu - FRANCE  
Tel. +33 (0)4 72 05 60 10 - Fax +33 (0)4 72 02 19 18 / [serf@serf.fr](mailto:serf@serf.fr)



**cv**.serf